

# Release Notes

## Analyst<sup>®</sup>

Software  
Version 1.5.2 HotFix 2

### Fixed Issue

#### **User-Created Formula columns in the Results Table are not automatically updated when samples are removed from the Results Table**

In the Quantitation module of the Analyst<sup>®</sup> 1.5.2 software, if one or more samples are removed from the Results Table, then user-created Formula columns in the Results Table are not automatically updated. As a result, the content in the Formula cells becomes out of sync with the sample entry rows that follow the deleted samples.

This incorrect data is presented in the Results Table. The data can be copied using Ctrl + C and then exported to a text or .pdf file, or printed. (AN-769)

### Where to Get Help

#### Support

- <http://sciex.com/request-support>

### Requirements

- Analyst<sup>®</sup> 1.5.2 with HotFixes to February 2011 Software

### Install the Software

This HotFix is compatible with Windows XP (SP 2 and SP3) operating system.

1. Log on to the computer as a user with Administrator privileges.
2. Deactivate the hardware profile and then close the Analyst software.
3. Extract the content of **Analyst152HF2.zip** file to your hard drive.
4. Navigate to this location on the hard disk, double-click the **setup.exe** file and then follow the on-screen instructions.

When the HotFix is installed, a new entry is added to the Programs and Features list in Control Panel. Additionally, a shortcut to these release notes is placed under **Start > All Programs > AB SCIEX > Analyst**.

#### Updated files

The HotFix adds the Analyst 1.5.2 HotFix 2 Release Notes.pdf to the:

- C:\Program Files\Analyst\Help folder.

The HotFix updates the QuantRT.ocx file located in the:

- C:\Program Files\Analyst\Bin folder

## Remove the HotFix

1. Deactivate the hardware profile and then close the Analyst<sup>®</sup> software.
2. Open **Control Panel** and then click **Add/Remove Programs**.
3. Select **Analyst<sup>®</sup> 1.5.2 HotFix 2** from the list and then click **Uninstall**.

The HotFix is removed from the program list. The Analyst 1.5.2 HotFix 2 Release Notes.pdf file and shortcut are also removed. The QuantRT.ocx file is reverted to its previous state.

## Known Issues and Limitations

The following patches are silently removed by the installation of **Analyst<sup>®</sup> 1.5.2 HotFix 2**:

- Analyst 1.5.2 Patch for Watson LIMS
- Analyst 1.5.2 Quant Patch for Cliquid

Contact support for update.

This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate those manufacturers' products as supplied by SCIEX for incorporation into its equipment and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are SCIEX's sole and exclusive representations, warranties, and obligations. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.

**For research use only. Not for use in diagnostic procedures.**

AB Sciex is doing business as SCIEX.

The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners.

AB SCIEX™ is being used under license.

© 2016 AB Sciex



AB Sciex Pte. Ltd.  
Blk 33, #04-06  
Marsiling Ind Estate Road 3  
Woodlands Central Indus. Estate  
SINGAPORE 739256

