

# LibraryView Software 1.8

## Release Notes

---



## Introduction

For information about a previous software version, refer to the document: *Release Notes* that came with that version of the software.

## Where to Get Help

To get access to all of the SCIEX software documentation, click **Start > SCIEX**.

### Software Guides and Tutorials

Document	Description
<i>LibraryView Software Quick Start Guide</i>	Supplies the procedures most frequently done by the user and a list of common terms and definitions.
<i>LibraryView Software Release Notes</i>	Supplies information about new software features, enhancements, and corrections, notes on use, known issues, and limitations.
<i>LibraryView Software Help</i>	Supplies the procedures for all of the tasks done by the user.

## New in Version 1.8

This section gives a description of the changes in the LibraryView software 1.8. To see the enhancements and corrected issues for an earlier version of the LibraryView software, refer to the document: *Release Notes* that came with that version of the software.

## New Features and Enhancements

Version 1.8 includes support for library searching of MS/MS spectra that were acquired in electron-activated dissociation (EAD) fragmentation mode:

- The properties for a transition for a compound include the fragmentation mode (**MSMS Fragmentation**) and, for transitions acquired with EAD, the kinetic energy value (**KE**). The value of **MSMS Fragmentation** is **CID** or **EAD**.
- The properties for an MS spectrum include **MSMS Fragmentation** and **KE**.
- The fragmentation mode is shown in the titles of MS/MS spectra. For EAD spectra, the kinetic energy (KE) is also included in the title.

- During the import of libraries created with earlier versions of the LibraryView software or the SCIEX OS software 3.4.5 or earlier, a fragmentation mode of collision-induced dissociation (CID) is automatically assigned to all of the MS/MS spectra.

## Known Issues and Limitations

- The new **MSMS Fragmentation** and **KE** fields are not exported to Excel format. (LBV-2330)

**Workaround:** Export to a different format.

- The new **MSMS Fragmentation** and **KE** fields are not printed on reports.
- Only the LibraryView software 1.6 and later versions can be upgraded to version 1.8. To upgrade earlier versions, contact SCIEX support.
- The SCIEX OS software 4.0 and LibraryView software 1.8 is the only supported combination if both are installed on the same computer.
- The installation of the LibraryView software does not complete if the user installs PostgreSQL, installs the LibraryView software, removes the LibraryView software, removes PostgreSQL, and then installs the LibraryView software again. (LBV-2148)
- During an export to an Excel file, if all of the columns are mapped to **---[not used]---**, then the software becomes unresponsive. (LBV-1561)
- If a spectrum that was acquired with an X500R QTOF system is added to the library, then the sample number in the graph title is incorrect. (LBV-1504)
- The Add User tab is not available in the standalone application. (LBV-217)
- The zooming feature has these issues (LBV-119):
  - Zooming is not available for graphs in the **List**, **Stack**, or **Grid** view.
  - If a graph is maximized, then it becomes unresponsive. To unlock the graph, click a sample name in the Acquired Spectrum pane.
  - If the Acquired Spectrum pane is maximized, then the zooming feature becomes unresponsive. To unlock the zooming feature, click in the Acquired Spectrum pane.
  - If the user changes between the TIC and XIC and then tries to change between Experiment 1 and Experiment 2 in the TIC pane, then the view might become locked and not change between the chromatograms in Experiment 1 and Experiment 2. To unlock the view, click the sample name in the Acquired Spectrum pane.
  - Spectral-like data might be shown instead of a chromatographic peak. To show the TIC correctly, click in the Acquired Spectrum pane to refresh the chromatogram.
- If a spectrum is added to a compound, then it is only shown in the **List** view. It is not shown in the **Grid** or **Stack** view. (LBV-116)

## Contact Us

### Addresses



AB Sciex Pte. Ltd.  
Blk33, #04-06 Marsiling Industrial Estate Road 3  
Woodlands Central Industrial Estate, Singapore 739256

#### SCIEX Headquarters

AB Sciex LLC  
250 Forest Street  
Marlborough, MA 01752  
USA

## Customer Training

- Global: [sciex.com/contact-us](https://sciex.com/contact-us)

## Online Learning Center

- [SCIEX Now Learning Hub](#)

## SCIEX Support

SCIEX and its representatives have a global staff of fully-trained service and technical specialists. They can supply answers to questions about the system or any technical issues that might occur. For more information, go to the SCIEX website at [sciex.com](https://sciex.com) or use one of the following links to contact us.

- [sciex.com/contact-us](https://sciex.com/contact-us)
- [sciex.com/request-support](https://sciex.com/request-support)

## Cybersecurity

For the latest guidance on cybersecurity for SCIEX products, visit [sciex.com/productsecurity](https://sciex.com/productsecurity).

## Documentation

This version of the document supersedes all of the previous versions of this document.

To see this document electronically, Adobe Acrobat Reader is required. To download the latest version, go to <https://get.adobe.com/reader>.

The latest versions of the documentation are available on the SCIEX website, at [sciex.com/customer-documents](https://sciex.com/customer-documents).

## LibraryView Software 1.8 Release Notes

---

---

**Note:** To request a free, printed version of this document, contact [sciex.com/contact-us](https://sciex.com/contact-us).

---

---

This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate such products as those manufacturers' products and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are the sole and exclusive representations, warranties, and obligations of SCIEX. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.

(GEN-IDV-09-10816-E)

For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd., or their respective owners, in the United States and/or certain other countries (see [sciex.com/trademarks](https://sciex.com/trademarks)).

AB Sciex™ is being used under license.

© 2025 DH Tech. Dev. Pte. Ltd.

