

DISCLAIMER

This software has not undergone extensive formal testing at this time. It will be included in the next version of the MPX™ Driver Software, at which time it will have been subjected to such testing. SCIEX assumes no liability or responsibility for any loss howsoever caused in connection with using this software.

MPX™ Driver

Software

Version 2.0

MPX™ Driver 2.0 Patch for No-Oven Configuration Support

This patch addresses the following issue:

In the MPX driver software version 2.0, the run will not start until the column oven reaches the temperature set in the method. This is to make sure that the run does not begin until the column oven is ready. This feature requires a column oven to be configured in the hardware profile.

This patch provides an option to run the MPX driver software samples as expected on a system without a column oven in the hardware configuration (Jira MPX -477).

If a column oven is present, then the oven temperature can be set to a value between 4 °C and 85 °C. If a column oven is not present, then the oven temperature must be set to 0 °C in the MPX driver software Method Editor.

Where to Get Help

Support

- <http://sciex.com/request-support>

Install the patch

1. Close the MPX™ driver software.
2. Deactivate the hardware profile and then close the Analyst software.
3. Double-click **MPX Driver 2.0 Patch for No-Oven Configuration Support.exe** file and then follow the instructions.

Updated files

The patch updates the following MPX™ driver software files:

C:\Program Files (x86)\SCIEX\MPX\Bin folder

- MPX.UI.exe

C:\Program Files(x86)\SCIEX\MPX\bin\Device Driver folder

- Multiplex Device Server.exe

Remove the patch

1. Open **Control Panel** and then click **Programs and Features**.
2. Select **MPX Driver 2.0 Patch for No-Oven Configuration Support** from the list and then click **Uninstall**.

Known Issues and Limitations

None

Notes on Use

None

This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate those manufacturers' products as supplied by SCIEX for incorporation into its equipment and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products, and are SCIEX's sole and exclusive representations, warranties, and obligations. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser, or for any adverse circumstances arising therefrom.

For research use only. Not for use in diagnostic procedures.

AB Sciex is doing business as SCIEX.

The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners.

AB SCIEX™ is being used under license.

© 2018 AB Sciex



AB Sciex Pte. Ltd.
Blk 33, #04-06
Marsiling Ind Estate Road 3
Woodlands Central Indus. Estate
SINGAPORE 739256

