

Molecule Profiler Software

1.3.4 Release Notes



Introduction

The Molecule Profiler software 1.3.4 is processing-only software, used for the identification and characterization of parent molecules and their metabolites or impurities.

The software processes accurate mass data, acquired from the following systems:

- TripleTOF systems, with the Analyst TF software 1.5 or later.
- X500 QTOF and ZenoTOF systems, with the SCIEX OS software.

Note: The numbers in parentheses are reference numbers for each issue or feature in the SCIEX internal tracking system.

New in Version 1.3.4

This section gives a description of the changes in the Molecule Profiler software 1.3.4. To see the enhancements and corrected issues for an earlier version of the Molecule Profiler software, refer to the document: *Release Notes* that came with that version of the software.

New Features and Enhancements

The Molecule Profiler software builds on the proven MetabolitePilot software. It includes these features and enhancements:

General

- The name of the **Custom Elements** feature was changed to **Chemical Dictionary**. To open the **Chemical Dictionary**, click **Setup > Chemical Dictionary**. Alternatively, the **Chemical Dictionary** can be opened directly from the Compound Library and Processing Parameters dialogs.
- ADC and peptide workflows: The algorithms that are used to assign the ion type for a fragment in which multiple isobaric species are possible have been improved.

Peptide workflow


- Added functionality to include data from CID and EAD fragmentation modes in one Results file.
- Added more functionality for the definition and analysis of cyclic peptides.
- Added more functionality to the **Chemical Dictionary** for the definition of neutral losses, radical losses, and unique ions in custom-defined amino acids.

Molecule Profiler Software 1.3.4 Release Notes

- Added a new biotransformation set, **Peptide Impurities**, for the identification of common impurities that can occur during peptide synthesis. This biotransformation set also includes the identification of these isobaric amino acid substitutions:
 - Leucine to isoleucine
 - Isoleucine to leucine
 - Aspartic acid to isoaspartic acid
 - Deamidation of asparagine to isoaspartic acid

Fixed Issues

General

Peptide and oligonucleotide workflows: In Interpretation view, for long sequences that extend across the whole pane, if the  **(Copy)** button or the right-click menu is used to copy the sequence, then some of the residues at the end of the sequence might not show when it is pasted. (MP-5500)

Peptide Workflow

On the Processing Parameters tab, if a peptide sequence that includes linkages and a building block with two or more capital letters in the name is identified in the Sequence pane, then the Catabolites table does not fill automatically. (MP-5497)

During results processing, if a peak can make a match to a y|a or y|b fragment, then the software might not automatically assign the correct fragment. (MP-5385)

During results processing, when the view mode is changed between Results view and Interpretation view, modified residues might not show in the correct color. (MP-5487)

During EAD results processing with processing parameters that contain only a CID reference spectrum, if an EAD spectrum is added with the paste function and then removed, then the calculation for the **Quality** score might not be correct. (MP-4345)

If a Results file is opened that contains an amino acid or amino acid modification that is not in the Chemical Dictionary, then the software might close unexpectedly. (MP-5486)

Structure and Sequence Assignment

- For lipid impurities that contain only an n-alkyl chain loss, for example a loss of n-butyl, the score is not calculated correctly for the proposed structure candidates. As a result, proposed candidates all have a rank of 1. Use the **Assigned Score** value as a guide for structure selection. (MP-5112)

Related Documentation

The documentation for the Molecule Profiler software is installed automatically with the software and is available from the Start menu: **SCIEX OS > SCIEX OS Documentation**.

Workflow-specific procedures are available through the **How Do I?** button in the top-right corner of the Molecule Profiler workspace in the SCIEX OS software. When **How Do I?** is clicked, users can select the Help topic from the list that is shown.

Installation

Notes on Installation

- The Molecule Profiler software 1.3.4 is installed as part of the SCIEX OS software, and activated with a version 1.3.1 license. For installation instructions and requirements for the SCIEX OS software, refer to the document: *SCIEX OS Software Installation Guide*. To activate the Molecule Profiler software, refer to the section: [Activate the Software](#).
- The Molecule Profiler software is removed with the SCIEX OS software. For instructions, refer to the document: *SCIEX OS Software Installation Guide*.
- To upgrade from an earlier version of the Molecule Profiler software, install the SCIEX OS software. The earlier version is removed and the new version is installed. Refer to the section: [Upgrade the Software](#).

Upgrade the Software

Note: In the peptide workflow, processing parameter files that were created in an earlier version of Molecule Profiler can be out of date after the upgrade. Open the parameter file to make sure that all of parameters are correct, and then save the file again before use.

Note: During installation, the processing parameter templates are overwritten with the newest versions of the templates. The user must reset any preferred default values and then save the default settings.

Note: Results files created in earlier versions of Molecule Profiler will be automatically recalculated when opened in this version and some values can change. The **MSMS Peak Area Assigned (%)** values in the data table can be different than the values that are displayed in the histograms. To update the histogram with the new values, process the data again.

Use this procedure to upgrade from earlier versions of the Molecule Profiler software to the Molecule Profiler software version 1.3.4.

Note: A Molecule Profiler software 1.3.1 license is required.

Molecule Profiler Software 1.3.4 Release Notes

1. Install the SCIEX OS software Refer to the document: *SCIEX OS Software Installation Guide*.
The installation program installs the SCIEX OS software and upgrades the Molecule Profiler software to version 1.3.4.
2. Activate the Molecule Profiler software version 1.3.4. Refer to the section: [Activate the Software](#).

Activate the Software

Note: To get the license, internet access is required. If the computer does not have Internet access, then make a copy of the generated computer ID. On a computer with Internet access, go to the licensing page of the SCIEX website, and then follow the instructions to get a license.

Note: Accept any changes from User Account Control during activation.

1. Open the SCIEX OS software.

Note: If the SCIEX OS software is not licensed, then the SCIEX OS Activation dialog opens. Go to step [4](#).

2. Open the Configuration workspace, and then go to the License page.
 3. Click **Install License**.
The SCIEX OS Activation dialog opens.
 4. Type the license key for the Molecule Profiler software in the applicable field.
The license key might be distributed on a printed activation certificate, or in an e-mail from SCIEX Now. If the license key is missing, then contact a SCIEX sales representative.
 5. Click **Copy ID to Clipboard**.
 6. Go to sciex.com/request-support.
 7. Follow the instructions to get the license.

After the required information is submitted, a license file is sent to all of the e-mail addresses supplied.
 8. Close the browser window.
 9. When the e-mail that contains the license file is received, copy the license file to the workstation desktop.
 10. In the SCIEX OS Activation dialog, click **Install License File**.
The Select the new license file to be installed dialog opens.
 11. Browse to and then select the license file.
 12. Click **Open**.
A confirmation dialog opens.
-

13. Click **OK**.

Note: Close the SCIEX OS software and then open it again. The Molecule Profiler tile is added to the Home page.

Known Issues and Limitations

General Issues

MetabolitePilot Software Compatibility

- In the Results and Correlation workspace, the average mass column for data files created in the MetabolitePilot 2.0.4 software contains **0** instead of **N/A**. (MP-2371)

Data Processing

- During data acquisition, if a user uses the Molecule Profiler software to process large amounts of data on the same computer, then the acquisition stops. To prevent this issue, do not process a large amount of data during data acquisition, or process the data on a separate computer.

Processing Method Issues

SWATH Acquisition MS/MS Reference Spectra

- During the extraction of reference spectra from SWATH acquisition data, the software proposes an extensive list of MS/MS spectral data. Some of the precursors related to the proposed MS/MS spectral data might have a low TOF MS peak intensity or a low chromatographic peak intensity. (MP-1854)

Product Ion and Neutral Losses Tab

Peptide and oligonucleotide workflows: If the user opens a processing method that does not have a spectrum on the Product Ion and Neutral Losses tab, adds a spectrum, and then clicks the **Assign Fragments** button, then the fragment table is not filled. To fill the table, change one of the filters, change it back, and then click **Assign Fragments**. (MP-3071)

Batch Workspace Issues

Peak Finding

- If more than one peak finding strategy is used to process a data file, then the chromatograms related to specific peak finding strategies might not be shown for some of the metabolites in the Results file. To make sure that all of the applicable chromatograms are shown, increase the **Maximum number of unexpected metabolites** on the MS Parameters tab of the Generic Parameters. (MP-2011)

- Peptide workflow, SWATH acquisition data: If an isotope pattern is used for peak finding, then only the singly-charged form of the fragment ion formula is used. (MP-2007)

Scheduled MRM (sMRM) Data Processing

- MS/MS data is not shown, and it is not used to calculate the score. (MP-2976)

Results Workspace Issues

Metabolite Name and Score

- For each metabolite, a list of possible MS identities is shown in the **Name** field of the Edit Name and Formula dialog. For ADC results, the MS identities that come from one or more antibody fragments with the same masses are not included in the list of other proposed names. Thus, the user cannot easily get access to them in the Interpretation view. (MP-1745)
- Different selections for **Source of Reference MS/MS Spectrum (Sample or Selected reference spectrum)** in the MS/MS Parameters of the Generic Parameters might show a different **MS/MS Similarity** score in the Details pane of the Results workspace. (MP-1839)

Grouping

- When the grouping feature is used, the header of the Results Table is not updated correctly after rows are deleted. As a workaround, click **Save**, remove all data from the group, and then add all data into the group again. The correct number is then shown. (MP-2929)

Interpretation

- For ADC results, the **Load Sequence** option automatically adds the protein fragment sequence that is related to the name assigned during data processing. This also true when the name of the metabolite has been changed with the **Edit Name and Formula** option. As a workaround, the sequence of interest can be typed in the Metabolite Sequence pane. (MP-1957)
- Oligonucleotide workflow: If the user applies a filter to the Fragments list that hides all of the fragments, and then clicks **Apply**, issues can occur. To prevent issues, make sure that the Fragments list contains at least one fragment, and then click **Apply**. (MP-3024)
- Oligonucleotide workflow: After a new MS/MS spectrum is added, and removal of interpretation data is confirmed, the **Assigned** check box is not cleared, and the sequence stays. The user can **Paste MS/MS** and **Assign Fragments**. (MP-3016)

Reports Issues

- Peptide or ADC workflows: Amino acid modifications that are present in assigned metabolite sequences and in sequences in the Fragments table are in square brackets in the software user interface. But, when an interpretation report is created for peptide and ADC results, the square brackets are not always included in the printed report. (MP-2186)

Contact Us

Addresses

SCIEX Headquarters	AB Sciex LLC 250 Forest Street Marlborough, MA 01752 USA
-------------------------------	---

Customer Training

- Global: sciex.com/contact-us

Online Learning Center

- [SCIEX Now Learning Hub](#)

SCIEX Support

SCIEX and its representatives have a global staff of fully-trained service and technical specialists. They can supply answers to questions about the system or any technical issues that might occur. For more information, go to the SCIEX website at sciex.com or use one of the following links to contact us.

- sciex.com/contact-us
- sciex.com/request-support

Cybersecurity

For the latest guidance on cybersecurity for SCIEX products, visit sciex.com/productsecurity.

Documentation

This version of the document supersedes all of the previous versions of this document.

The latest versions of the documentation are available on the SCIEX website, at sciex.com/customer-documents.

Note: To request a free, printed version of this document, contact sciex.com/contact-us.

This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate such products as those manufacturers' products and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are the sole and exclusive representations, warranties, and obligations of SCIEX. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.

(GEN-IDV-09-10816-E)

For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd., or their respective owners, in the United States and/or certain other countries (see sciex.com/trademarks).

AB Sciex™ is being used under license.

© 2025 DH Tech. Dev. Pte. Ltd.

