
DISCLAIMER

This software patch has not undergone extensive formal testing at this time. The patch has undergone sufficient testing to prove that it resolves the issue detailed below. This patch will be included in the next version of the ProteinPilot software, at which time it will have been subjected to complete and formal testing. SCIEX assumes no liability or responsibility for any loss howsoever caused in connection with using this software.

ProteinPilot 5.0.2 Patch for UniProt Issue (Rev.3)

This patch fixes the issue where UniProt annotations cannot be retrieved. This issue was corrected in version 5.0.2 of the ProteinPilot software and then in two subsequent patches. However, UniProt has made more changes on their website and the Uniprot annotations can no longer be retrieved using the ProteinPilot software (PP-828).

Install the Patch

Prerequisites
<ul style="list-style-type: none">The ProteinPilot software 5.0.2 is installed.

1. Log on to the computer as a user with Administrator privileges.
2. Wait for any acquisitions to complete.
3. Close the ProteinPilot software.
4. Double-click the `ProteinPilot-5.0.2-Patch-for-UniProt-Issue (Rev.3) .exe` file.
5. Follow the on-screen instructions to complete the installation.
Earlier revisions of this patch will be removed silently during the installation of this revision.

Updated Files

The ProteinPilot 5.0.2 Patch for UniProt Issue (Rev.3) replaces the following files in the `Program Files\AB SCIEX\ProteinPilot` folder:

- `progroup.exe`
- `ABSciex.ProteinObjects`

- `ABSciex.ProteinObjects_x64.dll`

Remove the Patch

1. Log on to the computer as a user with Administrator privileges.
2. Close the ProteinPilot software.
3. Open the Control Panel in Large icons or Small icons mode, and then click **Programs and Features**.
4. Select **ProteinPilot 5.0.2 Patch for UniProt Issue (Rev.3)** and then click **Uninstall**.

The patch is removed from the program list.

This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate those manufacturers' products as supplied by SCIEX for incorporation into its equipment and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are the sole and exclusive representations, warranties, and obligations of SCIEX. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.

(GEN-IDV-09-10816-D)

For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd., or their respective owners, in the United States and/or certain other countries (see sciex.com/trademarks).

AB Sciex™ is being used under license.

© 2023 DH Tech. Dev. Pte. Ltd.



AB Sciex Pte. Ltd.
Blk33, #04-06 Marsiling Industrial Estate Road 3
Woodlands Central Industrial Estate, Singapore 739256