
Natural Products HR-MS/MS Spectral Library

High Resolution Accurate Mass Libraries

Release Notes



This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate such products as those manufacturers' products and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are the sole and exclusive representations, warranties, and obligations of SCIEX. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.

(GEN-IDV-09-10816-E)

For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd., or their respective owners, in the United States and/or certain other countries (see sciex.com/trademarks).

AB Sciex™ is being used under license.

© 2024 DH Tech. Dev. Pte. Ltd.



AB Sciex Pte. Ltd.
Blk33, #04-06 Marsiling Industrial Estate Road 3
Woodlands Central Industrial Estate, Singapore 739256

Contents

1 Introduction	4
Features	4
Requirements	4
Supported Equipment	5
2 Notes on Use, Known Issues, and Limitations	6
3 Installation	7
Install a Licensed Natural Products HR-MS/MS Accurate Mass Library	7
Compound Conflicts	10
See Compound Conflicts	11
Merge Compounds	11
Overwrite Compounds	11
Keep Original Compounds	12
Update Libraries that Contain Added Spectra	12
Recommended Library Search Settings	12
Contact Us	14
Customer Training	14
Online Learning Center	14
SCIEX Support	14
Cybersecurity	14
Documentation	14

Features

The Natural Products HR-MS/MS Spectral Library 1.2 is an update to the Natural Products HR-MS/MS Spectral Library 1.1. This verified library includes High Resolution Accurate Mass Spectra for 1,295 natural products and traditional medicines. This library is released in the `lbp` format.

The library is compatible with:

- Data acquired with the SCIEX OS software 3.1 or later on a QTOF accurate mass system.

In the LibraryView software 1.4 or later and the SCIEX OS 3.1 software, use the library to do these tasks:

- Accurately identify compounds with more confidence in the reported results.
- Quickly search for compounds for both targeted and non-targeted screening.
- Use the comparative profiling option that is available in the SCIEX OS software to compare a sample to a control for qualitative review.

Requirements

SCIEX OS Software
An English (US) version of Windows 10 (64-bit) operating system.
The user must be logged on to the computer as a user with Administrator privileges.
Internet access is required to get a license file for each installed High Resolution Accurate Mass (HRAM) library.
At least one of the following applications is required to install a library: <ul style="list-style-type: none">• A licensed version of the SCIEX OS software, 3.1 or later.• A licensed version of the LibraryView (SCIEX OS) software, 1.4 or later.
A licensed version of the LibraryView (SCIEX OS) software, 1.4 or later, is required to edit the library.

Supported Equipment

- ZenoTOF 7600 accurate mass system
- X500 QTOF accurate mass system
- TripleTOF accurate mass system
- Dell Precision T3600 computer, or a later model, supplied by SCIEX, with a minimum of 32 GB of RAM

Notes on Use, Known Issues, and Limitations

2

Natural Products HR-MS/MS Spectral Library 1.2

The Natural Products HR-MS/MS Spectral Library 1.2 contains of 1,295 high-resolution spectral data on natural products and traditional medicines. It is an updated version of the Natural Products HR-MS/MS Spectral Library 1.1.

Sometimes the LibraryView software stops responding while it is loading the library, or results are not found during searches of the library

To correct these issues, do the following:

1. Make sure that the computer has 32 GB RAM available.
2. Start the computer again.
3. If the issue continues, then contact SCIEX Technical Support at sciex.com/request-support and request that the SQL databases be cleared and that the libraries be installed again.

Sometimes the SCIEX OS software stops responding when the Import compounds from LibraryView database feature is used

The Import compounds from LibraryView software database feature is used to import compound names and formulas from the library to the SCIEX OS software.

Install a Licensed Natural Products HR-MS/MS Accurate Mass Library

A licensed library can be installed from a `zip` file downloaded from the SCIEX website. The `zip` file contains compound names, compound transition information, and compound library spectra.

Note: Internet access is required to obtain the license.

Note: If the Natural Products HR-MS/MS Spectral Library 1.1 was installed, then download the Natural Products HR-MS/MS Spectral Library 1.2 and use the old license key.

1. Log on to the computer as a Windows user with administrator privileges.
2. Download the required `zip` file from the SCIEX website.

Tip! To prevent potential installation issues, save the file to a location other than the computer desktop.

3. After the download is complete, right-click the downloaded file, and then click **Extract All**.
4. Go to <https://sciex.com/support/activate-software>, and then log on with a SCIEX username and password.

Note: To create an account, follow the on-screen instructions to create an account.

The Activate Software page opens.

5. Select the applicable instrument in the **Select Your Instrument** field.

Tip! If the instrument is not listed, then go to the [SCIEX Now](#) profile for the logged-on user, and add the instrument information.

6. In the Windows search field, type `ipconfig /all` to get the physical (MAC) addresses of the computer.
A physical address, in the format `34-02-86-06-8A-05`, is shown for each active adapter.
7. Type all of the physical addresses in the **Computer ID** field.

Tip! A maximum of three physical addresses can be typed. Separate the addresses with a space. For example, `34-02-86-06-8A-05 34-02-86-06-8A-01 34-02-86-06-8A-09`.

Installation

8. In the the **License Key** field, type the license key. If a license key is not available, then contact sciex.com/contact-us.
The license key is distributed through an e-mail from [SCIEX Now](https://sciex.com). If the license key is not available, then contact a sciex.com/request-support sales representative.


Note: The key begins with the letters AID.

9. Click **Submit**.
After the required information is submitted, a license file is sent to the e-mail address registered to the [SCIEX Now](https://sciex.com) account.
10. Save the license file to the correct location:
- On a computer with version 1.4 of the LibraryView (SCIEX OS) software or later installed, save the license file in the C:/Program Files/SCIEX/LibraryView/LibraryViewFramework/Server folder.

After the license is added, the library importer dialog opens automatically.

Tip! As an alternative, to launch the library, run the `libraryviewpackager.exe` file located in the C:/Program Files/SCIEX/LibraryView/LibraryViewFramework/Packager folder.

11. Click the **LibraryView Package (*.lbp)** option on the Library Importer dialog.

Figure 3-1 Library Importer Dialog: Select Source Type for Import

12. Browse to the files extracted in step 3, and then select the TCM MSMS Library 2.2 lbp file.
13. To import all of the compounds, click **All** above the **Compound** column .

Figure 3-2 Library Importer Dialog: Select Compounds to Import

Library Name	Compound	Formula	CAS
TCM MS/MS Library 2.2	2,3,5,4-四羟基-2-苯基-2-呋喃葡萄糖 2,3,5,4'-Tetrahy...	C20H22O9	82373-94-2
	6-羟基-7,8-二甲氧基香豆素 6-hydroxy-7,8-dimet...	C11H10O5	
	(-)-丁香树脂酚-4-O-β-D-呋喃糖基- (1→2) -...	C33H44O17.HCOOH	136997-64-3
	(-)-丁香树脂酚-4-O-β-D-呋喃糖基- (1→2) -...	C33H44O17.Na	136997-64-3
	(-)-丁香树脂酚-4-O-β-D-呋喃糖基- (1→2) -...	C33H44O17	136997-64-3
	(-)-松脂素	C20H22O6	81446-29-9
	(+)-丁香脂素	C22H26O8	21453-69-0
	(+)-表松脂素-4'-O-β-D-葡萄糖苷	C26H32O11	24404-49-7
	(+)-表松脂素-4-O-β-D-葡萄糖苷	C26H32O11	74983-66-7
	(±)原苏木素B Protosappanin B	C16H16O6	102036-29-3
	(6aR, 11aR)-3-羟基-9, 10-二甲氧基紫檀烷 3-Hy...	C17H16O5	73340-41-7
	(7'R,8'S,8''R,8'''R)-表-丹酚酸B	C36H30O16	
	(R,S)-吉依春 Epigoltrin	C5H7NOS	1072-93-1
	1,2,3,6,7-五甲氧基山酮 1,2,3,6,7-Pentamethox...	C18H18O7	
	1,2,3,7-四甲氧基山酮 1,2,3,7-Tetramethoxyxan...	C17H16O6	22804-52-0
	1,4-二咖啡酰奎宁酸 1,4-Dicaffeoylquinic Acid (Cy...	C25H24O12	1182-34-9
	1,5-二咖啡酰奎宁酸 1,5-DCQA 1,5-Dicaffeoylqu...	C25H24O12	30964-13-7
	4,5-二甲氧基山酮 (4,5)-4,5-Dimethoxyxanth...	C15H12O4	

1295 from 1295 compounds selected.

- Click **Next**.

Note: If the import is cancelled before all of the compounds have been copied to the database, then any compounds that have already been imported stay in the database. The database will not go back to the pre-import state.

- For users with Natural Products HR-MS/MS Spectral Library 1.1, click **Overwrite All**. Refer to the section: [Compound Conflicts](#).
- Click **Finish**.
- If the LibraryView software contains any previous libraries, then use the **Export > Library as snapshot (lbp)** feature to make a backup of all of the libraries before installation.

Note: If issues occur with any of the libraries, then it is faster to import the entire library as a snapshot than to use the lbp packages to install all of the libraries again.

Compound Conflicts

When individual compounds or a library with a group of compounds is installed, the software searches the database for compounds with the same name or formula as each compound in the package. If compounds with the same name are found, then the software flags the compound in the package and shows a prompt to the user.

Users have the option to:

- Merge the compound information. New spectra, transitions, and retention times from the compound in the package are added to the compound information that is stored in the database.
- Overwrite the compound information. Compound information from the package replaces the compound information that is stored in the database.
- Keep compound information. Compound information in the database is kept and the compound information from the package is discarded.

Conflict information is available to help the user make the correct choice.

See Compound Conflicts

1. To see the details of the conflict, click **Resolve** beside the compound on the Library Importer dialog.
2. Do one of the following:
 - To keep the compound information that is in the database and discard the new information, click **Keep Original**.
 - To replace the compound information that is in the database with the new information, click **Use New**.
3. Do steps 1 and 2 for each compound.
4. After all of the conflicts are resolved, click **Finish**.

Merge Compounds

1. On the Library Importer dialog, do one of the following:
 - To merge new spectra, transitions, and retention times from individual compounds in the import package with the related compound information that is in the database, click **Merge**
 - To merge new spectra, transitions, and retention times from all of the compounds in the import package with the related compound information that is in the database, click **Merge All**.
2. After all of the conflicts are resolved, click **Finish**.

Overwrite Compounds

1. Do one of the following on the Library Importer dialog:
 - To overwrite all of the compound information that is in the database with the related compound information from the import package, click **Overwrite All**.

Installation

- To overwrite the compound information that is in the database with the related compound information from the import package, click **Resolve** beside the applicable compound and then click **Use New**.
2. After all of the conflicts are resolved, click **Finish**.

Keep Original Compounds

1. Do one of the following on the Library Importer dialog:
 - To keep all of the compound information that is in the database and discard the compound information from the import package, click **Keep All Original**.
 - To keep the individual compound information that is in the database and discard the compound information from the import package, click **Keep Original** beside the applicable compound.
2. After all of the conflicts are resolved, click **Finish**.

Update Libraries that Contain Added Spectra

1. To prevent the loss of customer-appended spectra when libraries are updated, do the following:
 - a. Append `compound_lab` name to the compound name.
 - b. Import the new library.
2. After all of the conflicts are resolved, click **Finish**.

Note: To delete compound spectra, a licensed version of LibraryView software is required.

Recommended Library Search Settings

The Natural Products HR-MS/MS Spectral Library 1.2 contains spectra from multiple instrument platforms and spectra acquired at different collision energy (CE) settings. Optimize the library search parameters carefully to effectively search the library. For the recommended library search settings for the SCIEX OS software, refer to the figure: [Figure 3-3](#). Complete the following settings:

- Deactivate the collision energy filtering to prevent false negatives.
- Set the **Precursor Mass Tolerance** and **Fragment Mass Tolerance** to 0.05 Da to decrease the processing time.
- For QTRAP system data, set **Precursor Mass Tolerance** and **Fragment Mass Tolerance** to 0.4 Da
- Reduce the **Intensity Threshold** to 0.02, if the fragment intensities fall below 0.05 (5%) of base peak

Figure 3-3 Recommended Search Settings for the Natural Products Library

Workflow
Components
Integration
Library Search
Calculated Columns
Flagging Rules
Advanced
Formula Finder
Non-targeted Peaks

Configure the library search parameters

☒ Perform Library Search

Library Search Algorithm: Candidate Search

Results Sorted By: Purity

Libraries To Search:
☒ Search All Libraries
☒ TCM MS/MS Library 2.2

Algorithm Parameters

☒ Precursor Mass Tolerance +/- 0.05 Da

☐ Collision Energy +/- 5 eV

☐ Retention Time +/- 0.5 min

Fragment Mass Tolerance +/- 0.05 Da

☐ Ignore Isotopes In Unknown

☐ Use Polarity

☐ Use Collision Energy Spread

☐ Use Compound Specific Purity Threshold

Maximal Number Of Hits: 5

Intensity Threshold: 0.01

Minimal Purity: 10.0 %

Intensity Factor: 5

Contact Us

Customer Training

- In North America: NA.CustomerTraining@sciex.com
- In Europe: Europe.CustomerTraining@sciex.com
- Outside the EU and North America, visit sciex.com/education for contact information.

Online Learning Center

- [SCIEX Now Learning Hub](#)

SCIEX Support

SCIEX and its representatives have a global staff of fully-trained service and technical specialists. They can supply answers to questions about the system or any technical issues that might occur. For more information, go to the SCIEX website at sciex.com or use one of the following links to contact us.

- sciex.com/contact-us
- sciex.com/request-support

Cybersecurity

For the latest guidance on cybersecurity for SCIEX products, visit sciex.com/productsecurity.

Documentation

This version of the document supersedes all previous versions of this document.

The latest versions of the documentation are available on the SCIEX website, at sciex.com/customer-documents.

Note: To request a free, printed version of this document, contact sciex.com/contact-us.
