
StatusScope[®] Remote Monitoring Service

Agent Upgrade Instructions



This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate those manufacturers' products as supplied by SCIEX for incorporation into its equipment and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are the sole and exclusive representations, warranties, and obligations of SCIEX. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.
(GEN-IDV-09-10816-B)

SCIEX Diagnostics products are for *in vitro* diagnostic use. Product(s) may not be available in all countries. For information on availability, please contact your local representative or refer to [sciex.com](https://www.sciex.com). All other SCIEX products are for research use only. Not for use in diagnostic procedures.

Trademarks and/or registered trademarks mentioned herein are the property of AB Sciex Pte. Ltd., or their respective owners, in the United States and/or certain other countries.

AB SCIEX™ is being used under license.

© 2020 DH Tech. Dev. Pte. Ltd.



AB Sciex Pte. Ltd.
Blk33, #04-06 Marsiling Industrial Estate Road 3
Woodlands Central Industrial Estate, Singapore 739256

Contents

1 Introduction.....	4
Features.....	4
Requirements.....	5
2 Add a User.....	6
Add a User.....	6
Add a User to the Administrator Group.....	9
Add a User to the Analyst® Software.....	13
Add a User to SCIEX OS	19
3 Installation.....	24
Prepare for Installation.....	24
Install the StatusScope® Remote Monitoring Service Agent.....	25
Troubleshooting Tips for Installation.....	31
Post-Installation Verification.....	32
Troubleshooting Tips for the Device Provisioning Task.....	35
4 Remove the StatusScope® Remote Monitoring Service Agent.....	39
A Supported Software.....	41
Analyst® Software Requirements.....	41
Analyst® TF Software Requirements.....	42
SCIEX OS Requirements.....	42
B Update Computers Configured with the Windows 7, 32-bit Operating System.....	43
C Enable the MSMQ Feature on Computers Configured with the Windows 7, 32-bit Operating System.....	44
Contact Us.....	46
Customer Training.....	46
Online Learning Center.....	46
SCIEX Support.....	46
CyberSecurity.....	46
Documentation.....	47

Introduction

1

This document is intended for StatusScope[®] Remote Monitoring Service customers who are upgrading an existing version of the software to the new version of the Agent. This update is required to enable customers to monitor their connected mass spectrometers through [SCIEX Now[™]](#) when we migrate to the new platform within the next few weeks.

The StatusScope[®] Remote Monitoring Service is used to remotely monitor SCIEX mass spectrometers to:

- Enhance system performance
- Enhance system serviceability
- Increase system uptime
- Enhance system use reporting capabilities
- Improve repair response time

Features

The StatusScope[®] Remote Monitoring Service uses a Web application infrastructure to receive performance parameters from, and report them to, connected mass spectrometers. With an established connection to a cloud server and a [SCIEX Now[™]](#) account, the status of the mass spectrometer can be monitored.

The StatusScope[®] Remote Monitoring Service provides:

- Support for multiple types of mass spectrometers and LC systems.
- A unified view of important system performance parameters for monitored mass spectrometers and their connected peripherals.
- Monitoring of the sample queue for jobs submitted through the Analyst[®] Software, the Analyst[®] TF Software, and SCIEX OS, on all supported mass spectrometers.
- A single access point for customers, through [SCIEX Now[™]](#), where customers can monitor the status, alerts, and notifications for all mass spectrometers that are connected to the StatusScope[®] Remote Monitoring Service. [SCIEX Now[™]](#) is configured for on-the-go viewing using a mobile device.
- Custom reports on mass spectrometer utilization, inventory, and service history.

Requirements

The new version of the StatusScope[®] Remote Monitoring Service Agent requires:

- A dedicated internet connection for the computer connected to the mass spectrometer, so that the StatusScope[®] Remote Monitoring Service can function as intended.
- An active [SCIEX Now[™]](#) account.
- A minimum of 2 GB free space on the C:/ drive.
- Outbound connectivity on HTTPS and MQTT protocols, port 443. The firewall must allow outbound communication from StatusScope[®] Remote Monitoring Service to aws-iot, including:
 - a191fe2mx36xf8-ats.iot.us-east-1.amazonaws.com
 - iotregistrationcertificates.s3.amazonaws.com

Contact a local IT network administrator to confirm that this requirement is met.

Add a User

2


The procedures in this chapter must be performed on the acquisition computer. A user who is part of the Administrator group for Windows, and who is assigned to the Administrator role in the Analyst® Software and SCIEX OS, if applicable, must complete the installation.

Note: All of the dialogs provided in this chapter were captured on a computer configured with the Windows 10 operating system. Comparable dialogs are shown on computers configured with the Windows 7, 32-bit or 64-bit operating system. The dialogs are provided as examples only.

Add a User


1. Right-click **Start > Search** and then type **computer management** in the **Search** field.

Figure 2-1 Computer Management App




2. Click the **Computer Management App**.

Figure 2-2 Computer Management Window



3. Click **Local Users and Groups**.

Figure 2-3 Local Users and Groups Options



4. Right-click **Users** and then click **New User**.
The New User dialog opens.
5. Complete the dialog as follows:
 - a. Clear the **User must change password at next logon** check box.


Add a User

- b. Type a **User name** in the field provided. For example, **StatusScopeUser**.

Note: The user name cannot contain spaces.

- c. Type a **Password** in the field provided.
- d. Type the same password in the **Confirm Password** field.

Figure 2-4 New User Dialog




The 'New User' dialog box is shown with the following fields and options:

- User name:** StatusScopeUser
- Full name:** (empty)
- Description:** (empty)
- Password:** (masked with dots)
- Confirm password:** (masked with dots)
- ☐ User must change password at next logon
- ☐ User cannot change password
- ☐ Password never expires
- ☐ Account is disabled
- Buttons:** Help, Create, Close

6. Click **Create**.
7. Click **Close**.
8. On the Computer Management window, double-click **Users**.
The name of the new user is shown in the list of users.

Figure 2-5 Computer Management Window – List of Users




9. Continue with [Add a User to the Administrator Group](#) on page 9.

Add a User to the Administrator Group

1. On the Computer Management window, double-click **Groups**.


Figure 2-6 Computer Management Groups Options Window



Add a User


2. Right-click **Administrators** and then click **Add to Group**.

Figure 2-7 Administrators Properties Dialog



3. Click **Add**.


Figure 2-8 Select Users, Computers, Service Accounts, or Groups Dialog



4. Make sure that the correct location is shown in the **From this location** field.

Note: The location is typically the computer name. To change the selection, click **Locations**, select the appropriate location on the Locations dialog, and then click **OK**.


Figure 2-9 Locations Dialog



5. Type the name of the user that was created in [Add a User on page 6](#) in the **Enter the object names to select** field and then click **Check Names**. For example, **StatusScopeUser**.

The name is shown, underlined, in the **Enter the object names to select** field.


Figure 2-10 Select Users Dialog – Valid User



Add a User


Tip! If the name of the user cannot be verified, then make sure that the name has been typed correctly and that the user is a valid computer or network domain user.

Figure 2-11 Name Not Found Dialog



-
6. Click **OK** to close the Select Users dialog.

The name of the new user is shown in the **Members** field.

Figure 2-12 Administrators Properties Dialog – Administrator Added

7. Click **Apply** and then click **OK**.
8. Close the Computer Management dialog.
9. Add the user to the Administrator role in the acquisition software:
 - [Add a User to the Analyst® Software on page 13](#)
 - [Add a User to SCIEX OS on page 19](#)

Add a User to the Analyst® Software

Note: In this procedure, all references to the Analyst® Software can be replaced with the Analyst® TF Software.


Note: Any changes to the Analyst® Software security configuration take effect after the Analyst® Software is restarted.

1. Open the Analyst® Software.

Add a User

2. On the Navigation bar, under **Configure**, double-click **Security Configuration**.


Figure 2-13 Security Configuration Dialog



3. Click the **People** tab.


Note: If the People tab is not shown, then the **Security Mode** in the Analyst[®] Software is set to **Single User Mode** and this procedure does not have to be completed. Continue with [Prepare for Installation on page 24](#).

Figure 2-14 Security Configuration Dialog – People Tab



4. Click **New Person**.

Figure 2-15 Select Users or Groups Dialog




5. Make sure that the correct location is shown in the **From this location** field.

Note: The location is typically the computer name. To change the selection, click **Locations**, select the appropriate location on the Locations dialog, and then click **OK**.

Add a User


Figure 2-16 Locations Dialog



6. Type the name of the user that was created in [Add a User on page 6](#) in the **Enter the object names to select** field and then click **Check Names**. For example, **StatusScopeUser**.


The name is shown, underlined, in the **Enter the object names to select** field.

Figure 2-17 Select Users or Groups Dialog - User Added




Tip! If the name of the user cannot be verified, then make sure that the name has been typed correctly and that the user is a valid computer or network domain user.

Figure 2-18 Name Not Found Dialog



7. Click **OK**.

Figure 2-19 Security Configuration Dialog




Add a User

- Make sure that the name of the user that was added is selected and then click **Administrator** in the **Available Roles** field and then click .


Administrator is shown in the **Role(s) Selected** field.

Figure 2-20 Security Configuration Dialog



- Click **Apply**.

Figure 2-21 SecurityConfig Message



- Click **OK** to close the message.
- Click **OK** to close the Security Configuration dialog.
- Close the Analyst® Software.


Note: Any changes to the Analyst® Software security configuration take effect after the Analyst® Software is restarted.

13. Continue with [Prepare for Installation on page 24](#).

Add a User to SCIEX OS

1. Open SCIEX OS.


Figure 2-22 SCIEX OS Home Page



2. Open the Configuration workspace.


Add a User

Figure 2-23 Configuration Workspace




3. Click **User Management**.

Figure 2-24 User Management




4. Click **Add User** ().

Figure 2-25 Select User or Group Dialog


5. Make sure that the correct location is shown in the **From this location** field.

Note: The location is typically the computer name. To change the selection, click **Locations**, select the appropriate location on the Locations dialog, and then click **OK**.

Figure 2-26 Locations Dialog

6. Type the name of the user that was created in [Add a User on page 6](#) in the **Enter the object names to select** field and then click **Check Names**. For example, **StatusScopeUser**.


The name is shown, underlined, in the **Enter the object names to select** field.

Figure 2-27 Select User or Group Dialog – User Added

Add a User


Tip! If the name of the user cannot be verified, then make sure that the name has been typed correctly and that the user is a valid computer or network domain user.

Figure 2-28 Name Not Found Dialog




7. Click **OK**.

Figure 2-29 User Management – User Added



8. Make sure that the name of the user that was added is shown in the User name or group field and then select the **Active user or group** check box.
9. In the **Roles** area, select the check box to the left of **Administrator**.

Figure 2-30 User Management – Administrator Selected



10. Click **Save**.
11. Close SCIEX OS.
12. Continue with [Prepare for Installation on page 24](#).

Prepare for Installation

Prerequisites

- Make sure that the user who will be completing the installation is assigned to the Administrator group. Refer to [Add a User to the Administrator Group on page 9](#).
- Make sure that the user who will be completing the installation is assigned to the Administrator role in the acquisition software. Refer to [Add a User to the Analyst® Software on page 13](#) or [Add a User to SCIEX OS on page 19](#).
- (For computers configured with the Windows 7, 32-bit operating system only) Make sure that SP1 and the Security Update for Windows 7 (KB2813430) are installed. Refer to [Update Computers Configured with the Windows 7, 32-bit Operating System on page 43](#).

1. Download the required zip file from sciex.com/software-support/software-downloads.
2. Save the file to a local drive on the computer that is connected to the instrument to be monitored.

Tip! To prevent potential installation issues, save the file to a location other than the computer desktop and disconnect any external USB storage devices.

3. After the download is complete, right-click the downloaded file and then click **Extract All** to extract the installation package.

Note: By default, the files are extracted to the location where the zip file was saved. Make sure that the path to the installation folder is not too long. If the path is longer than 118 characters, then installation might not proceed.

4. Continue with [Install the StatusScope® Remote Monitoring Service Agent on page 25](#).

Install the StatusScope® Remote Monitoring Service Agent

Prerequisites


- The installation must be performed on the acquisition computer.

IMPORTANT:

If a previous version of the StatusScope® Remote Monitoring Service is currently installed on the computer, it should not be removed. The StatusScope® Remote Monitoring Service is required to enable users to continue monitoring the status of their systems. Consult with a Field Service Employee (FSE) before removing the existing software.

1. Log on to the computer as the user that was created in [Add a User on page 6](#).
2. Make sure that the acquisition software is open.
3. Make sure that the hardware profile in the acquisition software is active and that it contains the mass spectrometer that is being monitored.
4. Navigate to the installation folder, **StatusScopeSetup-2.0**, open the folder, and then open the **Install > Utilities** folder.
5. Double-click **StatusScopePortCheck.exe**.

Figure 3-1 Connectivity Confirmation



6. Make sure that **Yes** is shown in the **Open** column for both servers.


Installation

Tip! Contact a local IT network administrator to confirm that the firewall is configured to allow outbound communication from the StatusScope® Remote Monitoring Service to the two host sites shown in [Figure 3-1](#). Refer to [Requirements on page 5](#).


7. Press **Enter** to close the dialog.
 8. Return to the root **StatusScopeSetup-2.0** folder, right-click **StatusScopeSetup.exe** and then select **Run as administrator**.
-

Tip! If the computer is configured with version 1607 or 1703 of the Windows 10 operating system, then the following message is shown. Continue with step [9](#).


Figure 3-2 .Net Framework Installation Dialog



Tip! If the computer is configured with any other version of the Windows 10 operating system, then the following dialog is shown. Continue with step [12](#).


Figure 3-3 End User License Agreement

Tip! If the computer is configured with the Windows 7, 32-bit operating system, then the following message is shown. Continue with [Enable the MSMQ Feature on Computers Configured with the Windows 7, 32-bit Operating System on page 44](#).

Figure 3-4 Install MSMQ Message Dialog


9. Click **OK**.

Figure 3-5 .NET Framework Installation Complete Dialog



10. Click **Finish**.

Figure 3-6 .NET Framework Restart Computer Dialog




11. Click **Restart Now**.

Note: When the computer has started again, return to step [1](#) of this procedure.

12. Click **I Agree**.


Note: It might take a few minutes for the next screen in the sequence to be shown. No user intervention is required.

Figure 3-7 Administrator User Credentials

The screenshot shows a window titled "SCIEX StatusScope Agent" with a close button (X) in the top right corner. The main text reads: "Please provide the user details to run the Statusscope Agent". Below this, there are three bullet points: "• User should have Administrative rights.", "• User should have permission to run Windows services.", and "• User should be configured in the Acquisition software with necessary permissions." Underneath the bullet points, there are two input fields: "Username:" and "Password:". At the bottom of the window, there are three buttons: "Previous", "Next", and "Cancel".

13. Type the name of the user that was created in [Add a User on page 6](#) in the **Username** field.
14. Type the **Password** in the field provided and then click **Next**.

Note: It might take a few minutes for the next screen in the sequence to be shown. No user intervention is required.


Figure 3-8 Installation Progress

The screenshot shows a window titled "SCIEX StatusScope Agent" with a close button (X) in the top right corner. The main text reads: "Current component: Microsoft.VisualStudio.C++.14.15.x64, Version : 14.0.23026.0". Below this, it says "Overall Progress :". Underneath, there is a progress bar that is partially filled with green. At the bottom right of the window, there is a "Cancel" button.

15. If both the Analyst[®] Software or Analyst[®] TF Software and SCIEX OS are installed on the workstation, then select the **Agent Type** that matches the software that is used for acquisition.

Note: If only one of the software applications is installed, then this dialog is not shown. Both Analyst[®] and Analyst[®] TF Software are shown as Analyst in the **Agent Type** field.


Figure 3-9 Agent Type



16. Click **Next**.

Note: It might take a few minutes for the next screen in the sequence to be shown. No user intervention is required.

Figure 3-10 Installation Complete





17. Click **Finish**.

Note: The system might take a few minutes to initialize and for the dialog to close. No user intervention is required.

18. Continue with [Post-Installation Verification on page 32](#).



Troubleshooting Tips for Installation

Table 3-1 Installation Issues

Issue	Resolution
	Install one of the following acquisition software: <ul style="list-style-type: none">Analyst® Software. Refer to Analyst® Software Requirements on page 41 for a list of supported versions.Analyst® TF Software. Refer to Analyst® TF Software Requirements on page 42 for a list of supported versions.SCIEX OS. Refer to SCIEX OS Requirements on page 42 for a list of supported versions.
	Make sure that the user is assigned to the Administrator group. The user account can be local to the computer or part of a network domain. Refer to Add a User to the Administrator Group on page 9 .

Installation

Table 3-1 Installation Issues (continued)


Issue	Resolution
	Make sure that the username and the password are typed correctly.
	Send an e-mail to StatusScope.Migration@sciex.com . The e-mail must include the contact information and the serial number of the connected mass spectrometer.

Post-Installation Verification

Note: All of the dialogs provided in this procedure were captured on a computer configured with the Windows 10 operating system. Comparable dialogs are shown on computers configured with the Windows 7, 32-bit or 64-bit operating system. The dialogs are provided as examples only.


1. Right-click **Start > Search** and then type **services** in the **Search** field.

Figure 3-11 Services App



2. Click the **Services App**.

Figure 3-12 Services Window



3. Verify that the status of the **StatusScopeAnalystAgent** or the **StatusScopeSciexOSAgent** service, and the **StatusScopeFieldGateway** service is **Running**.

Note: On computers configured with the Windows 7 operating system, the status of the service is **Started**.

Installation

Figure 3-13 StatusScope Services

StatusScope Access Remote	Provides se...	Manual
StatusScopeAnalystAgent	StatusScope...	Running Automatic
StatusScopeFieldGateway	Field gatew...	Running Automatic
StatusScopeSciexOSAgent	StatusScope...	Automatic

4. Verify that the installation logs are available in the C:\Users\[username]\AppData\Local\Temp folder, where [username] is the name of the user that is logged on the computer and has completed the installation.

Figure 3-14 Installation Logs

Setup_20200505074809.log	2020/05/05 7:55 AM	Text Document	1 KB
SCIEX_StatusScope_Agent_20200505074812.log	2020/05/05 7:55 AM	Text Document	460 KB
PDApp.log	2020/05/05 7:55 AM	Text Document	9,591 KB
SCIEX_StatusScope_Agent_20200505074812_007_NotifierSetup.log	2020/05/05 7:54 AM	Text Document	129 KB
SCIEX_StatusScope_Agent_20200505074812_006_FieldGatewaySetup.log	2020/05/05 7:54 AM	Text Document	109 KB
SCIEX_StatusScope_Agent_20200505074812_005_SciexOSAgentSetup.log	2020/05/05 7:54 AM	Text Document	365 KB
SCIEX_StatusScope_Agent_20200505074812_004_AnalystAgentSetup.log	2020/05/05 7:54 AM	Text Document	148 KB
18e190413af045db88dfbd29609eb877.db.ses	2020/05/05 7:54 AM	SES File	1 KB
dd_vcrist_x86_20200505075403.log	2020/05/05 7:54 AM	Text Document	8 KB
mat-debug-30072.log	2020/05/05 7:54 AM	Text Document	0 KB
dd_vcrist_amd64_20200505075401.log	2020/05/05 7:54 AM	Text Document	7 KB

5. Verify that the application logs are available in the C:\ProgramData\SCIEX\Logs\StatusScope folder.

Figure 3-15 Application Logs

Name	Date modified	Type	Size
StatusScopeDeviceAgent.log	2020/05/05 8:00 AM	Text Document	99 KB
StatusScopeFieldGateway.log	2020/05/05 8:04 AM	Text Document	355 KB
StatusScopeNotifier.log	2020/05/05 7:54 AM	Text Document	3 KB

6. Do one of the following:
 - On computers configured with the Windows 10 operating system, navigate to the C:\Program Files (x86)\SCIEX\StatusScope\Utilities folder, right-click **StatusScopeTest.exe**, and then select **Run as administrator**.
 - On computers configured with the Windows 7 operating system, navigate to the C:\Program Files\SCIEX\StatusScope\Utilities folder, right-click **StatusScopeTest.exe**, and then select **Run as administrator**.

This test tool validates that all prerequisites, services, data in message queue, and provisioning has been completed.

Figure 3-16 StatusScope Agent Troubleshooter Window – Device Provisioning Successful

```

=====
StatusScopeFieldGateway service is installed: True and it is in Running state : True
=====
Agent Queue contains 1 messages
Normal Queue contains 0 messages
=====
Java path is set in environment PATH variable.
Java path is: C:\Program Files\Amazon Corretto\jdk1.8.0_242\bin
Mosquitto path is set in environment PATH variable.
Mosquitto path is: C:\Program Files\mosquitto
OpenSSL 64-bit path is set in environment PATH variable.
OpenSSL 64-bit path is: C:\Program Files\OpenSSL-Win64\bin
cURL path is set in environment PATH variable.
cURL path is: C:\Program Files\cURL\bin
=====
Checking for Device Provisioning....
=====
Device Provisioning is Successful.
Thing name is: SIM-f103f82e-c34f-47cc-b78b-4edce8f1e5b7
Thing group name is: QTRAP4500
=====

```

7. Make sure that the statement **Device Provisioning is Successful** is shown at the bottom of the StatusScope Agent Troubleshooter window. If the statement **Device Provisioning Failed** is shown, then continue with [Troubleshooting Tips for the Device Provisioning Task on page 35](#).

Figure 3-17 StatusScope Agent Troubleshooter Window – Device Provisioning Failed

```

=====
StatusScopeFieldGateway service is installed: True and it is in Running state : True
=====
Agent Queue contains 1 messages
Normal Queue contains 0 messages
=====
Java path is set in environment PATH variable.
Java path is: C:\Program Files\Amazon Corretto\jdk1.8.0_242\bin
Mosquitto path is set in environment PATH variable.
Mosquitto path is: C:\Program Files\mosquitto
OpenSSL 64-bit path is set in environment PATH variable.
OpenSSL 64-bit path is: C:\Program Files\OpenSSL-Win64\bin
cURL path is set in environment PATH variable.
cURL path is: C:\Program Files\cURL\bin
=====
Checking for Device Provisioning....
=====
Thing name is:
Thing group name is:
Device Provisioning failed.
=====


```

8. Press **Enter** to exit the dialog.

Troubleshooting Tips for the Device Provisioning Task


1. Right-click **Start > Search** and then type **services** in the **Search** field.

Figure 3-18 Services App



2. Click the **Services App**.

Figure 3-19 Services Window



3. Verify that the status of the **StatusScopeAnalystAgent** or the **StatusScopeSciexOSAgent** service, and the **StatusScopeFieldGateway** service is **Running**.

Note: On computers configured with the Windows 7 operating system, the status of the service is **Started**.

Figure 3-20 StatusScope Services

StatusScope Access Remote	Provides se...	Manual
StatusScopeAnalystAgent	StatusScope...	Running Automatic
StatusScopeFieldGateway	Field gatew...	Running Automatic
StatusScopeSciexOSAgent	StatusScope...	Automatic

- For each service, right-click the service name and then click **Restart**.

Note: If one or both of the services are not running and cannot be started manually, then restart the computer and return to [Post-Installation Verification on page 32](#).

- Verify that the **Log On As** field for the **StatusScopeAnalystAgent** service or the **StatusScopeOSAgent** service contains the name of the user that was created in [Add a User on page 6](#). For example, **StatusScopeUser**.

Figure 3-21 StatusScopeAnalystAgent Log On As User

Name	Description	Status	Startup Type	Log On As
StatusScopeAnalystAgent	StatusScope agent service for Anal...	Running	Automatic	.\StatusScopeUser
StatusScopeFieldGateway	Field gateway service for commun...	Running	Automatic	Local System

- Do one of the following:
 - On computers configured with the Windows 10 operating system, navigate to the C:\Program Files (x86)\SCIEX\StatusScope\FieldGateway\JITP\JITP_Device_Register folder.
 - On computers configured with the Windows 7 operating system, navigate to the C:\Program Files\SCIEX\StatusScope\FieldGateway\JITP\JITP_Device_Register folder.

Figure 3-22 JITP_Device_Register Folder File List

Name	Date modified	Type	Size
chk_hst.txt	2020/06/04 5:46 PM	Text Document	5 KB
deviceCert.crt	2020/06/04 5:46 PM	Security Certificate	2 KB
deviceCert.csr	2020/06/04 5:46 PM	CSR File	2 KB
deviceCert.key	2020/06/04 5:46 PM	KEY File	2 KB
deviceCertAndCACert.crt	2020/06/04 5:46 PM	Security Certificate	3 KB
initiate_device_registration.bat	2020/05/25 11:37 PM	Windows Batch File	5 KB
properties.yaml	2020/06/04 5:46 PM	YAML File	1 KB
root.cert	2020/06/04 5:46 PM	CERT File	2 KB
rootCA.key	2020/06/04 5:46 PM	KEY File	2 KB
rootCA.pem	2020/06/04 5:46 PM	PEM File	2 KB
rootCA.srl	2020/06/04 5:46 PM	SRL File	1 KB

Installation

Note: If the **JITP_Device_Register** folder contains all of the files shown in [Figure 3-22](#), then the provisioning was successful. If the **JITP_Device_Register** folder contains more files or fewer files than those shown, then the provisioning did not occur or it failed.


7. Navigate to C:\ProgramData\SCIEX\Logs>StatusScope and download the **DeviceprovisioningLog.txt** file.
8. Send an email containing the following information to StatusScope.Migration@sciex.com:
 - contact information
 - serial number of the connected mass spectrometer
 - **DeviceprovisioningLog.txt** file

Remove the StatusScope[®] Remote Monitoring Service Agent

4


1. Log on to the computer as a user with Administrator privileges.
2. Click **Start > Control Panel**.
3. In the **Large icons** or **Small icons** view, click **Programs and Features**.
4. Right-click **SCIEX StatusScope Agent** and then click **Uninstall**.

Figure 4-1 Removal Wizard



5. Follow the on-screen instructions.

Figure 4-2 Removal Complete



6. Click **Finish**.

Analyst[®] Software Requirements

The following versions of the Analyst[®] Software are supported by the StatusScope[®] Remote Monitoring Service, including all of the available add-ons (Patches, HotFixes and Components releases):

- Version 1.6
- Version 1.6 Software Components for 6500 Series Instruments
- Version 1.6.1
- Version 1.6.2
- Version 1.6.2 HotFix to February 2015
- Version 1.6.2 Software Components for 3500 Series Instruments
- Version 1.6.3
- Version 1.6.3 HotFix 1
- Version 1.6.3 HotFix 2
- Version 1.6.3 HotFix 3
- Version 1.6.3 HotFix 4
- Version 1.7
- Version 1.7 with HotFix 1
- Version 1.7 HotFix 2
- Version 1.7 HotFix 3
- Version 1.7.1
- Version 1.7.1 HotFix 1

Analyst[®] TF Software Requirements

The following versions of the Analyst[®] TF Software are supported by the StatusScope[®] Remote Monitoring Service, including all of the available add-ons (Patches, HotFixes and Components releases):

- Version 1.7
- Version 1.7.1
- Version 1.7.1 HotFix 1
- Version 1.7.1 With Components for LC Devices
- Version 1.7.1 With Components for NanoCell
- Version 1.8
- Version 1.8 HotFix 1
- Version 1.8.1

SCIEX OS Requirements

The following versions of SCIEX OS are supported by the StatusScope[®] Remote Monitoring Service, including all of the available add-ons (Patches, HotFixes and Components releases):

- Version 1.6.1
- Version 1.6.2
- Version 1.7

Update Computers Configured with the Windows 7, 32-bit Operating System

B


1. Download SP1 for the Windows 7, 32-bit operating system from <https://www.microsoft.com/en-us/download/confirmation.aspx?id=5842>.
2. Follow the instructions to install SP1.
3. Download the Security Update for Windows 7 (KB2813430) from <https://www.microsoft.com/en-us/download/details.aspx?id=39110>.
4. Follow the instructions to install the security update.

Enable the MSMQ Feature on Computers Configured with the Windows 7, 32-bit Operating System

C


If the computer where the StatusScope® Remote Monitoring System Agent is being installed is configured with the Windows 7, 32-bit operating system, then the following message is shown when the **StatusScopeSetup.exe** file is double-clicked.

Figure C-1 Install MSMQ Message Dialog



1. Click **OK** to close the message.
2. Click **Start > Control Panel**.
3. In the **Large icons** or **Small icons** view, click **Programs and Features**.
4. In the left navigation panel, click **Turn Windows features on or off**.


Figure C-2 Programs and Features Window



Enable the MSMQ Feature on Computers Configured with the Windows 7, 32-bit Operating System

5. Select the **Microsoft Message Queue (MSMQ) Server** check box.

Figure C-3 Windows Features Dialog



6. Click **OK**.
7. Close the Programs and Features dialog.
8. Restart the computer and then continue with [Install the StatusScope® Remote Monitoring Service Agent on page 25](#).

Contact Us

Customer Training

- In North America: NA.CustomerTraining@sciex.com
- In Europe: Europe.CustomerTraining@sciex.com
- Outside the EU and North America, visit sciex.com/education for contact information.

Online Learning Center

- [SCIEX University™](#)

SCIEX Support

SCIEX and its representatives maintain a staff of fully-trained service and technical specialists located throughout the world. They can answer questions about the system or any technical issues that might arise. For more information, visit the SCIEX website at sciex.com or contact us in one of the following ways:

- sciex.com/contact-us
- sciex.com/request-support

Note: For any questions or issues related to the StatusScope® Remote Monitoring Service, select **Instrument & Hardware** as the **Issue Type**.

- StatusScope.Migration@sciex.com

CyberSecurity

For the latest guidance on cybersecurity for SCIEX products, visit sciex.com/productsecurity.

Documentation

This version of the document supercedes all previous versions of this document.

To view this document electronically, Adobe Acrobat Reader is required. To download the latest version, go to <https://get.adobe.com/reader>.

The latest versions of the documentation are available on the SCiEX website, at sciex.com/customer-documents.

Note: To request a free, printed version of this document, contact sciex.com/contact-us.
