AB SCIEX Reporter Template Suite 3.2 Component for Office 2024 Release Notes

Introduction

For information about an earlier software version, refer to the document: *Release Notes* that came with that software version.

Where to Get Help

- *Reporter Software 3.2 User Manual*
- *Reporter Software 3.2 Release Notes*

New in AB SCIEX Reporter Template Suite 3.2 Component for Office 2024

Templates that are used by the Reporter software can now be edited in Microsoft Office 2024. AB SCIEX Reporter Template Suite 3.2 Component for Office 2024 can be installed only if AB SCIEX Reporter Template Suite 3.2 is installed.

Known Issues and Limitations

After the AB SCIEX Reporter Template Suite 3.2 Component for Office 2024 is installed and before a template is edited, use the Reporter Template Updater utility in the `C:\Program Files\AB SCIEX\Reporter Template Updater` folder to update the template for the Reporter software to a compatible format. If users try to edit templates before the templates are updated with the Reporter Template Updater utility, then an error message is shown.

Installation

Install the Component for Office 2024

Prerequisites
• AB SCIEX Reporter Template Suite 3.2 is installed.

Note: AB SCIEX Reporter Template Suite 3.2 Component for Office 2024 is installed automatically during the installation of other SCIEX software products that include AB SCIEX Reporter Template Suite 3.2 Component for Office 2024, for example, the SCIEX OS software 5.0.1.

Note: If AB SCIEX Reporter Template Suite 3.2 HotFix 1 is installed, then it will be removed silently during the installation of AB SCIEX Reporter Template Suite 3.2 Component for Office 2024.

1. To let users edit templates in the Reporter software, install the Template Content Control Manager utility:
   a. Go to the `C:\Program Files\AB SCIEX\ReporterOfficeAddins\TemplateContentControlManager` folder.
   b. Double-click the `TemplateContentControlManager.vsto` file.
2. Update the templates to a compatible format. Refer to the section: Update the Templates for the Reporter Software.

Update the Templates for the Reporter Software

Use the Reporter Template Updater utility to update templates for the Reporter software to a format that is compatible with AB SCIEX Reporter Template Suite 3.2 Component for Office 2024.

1. Go to the `C:\Program Files\AB SCIEX\ReporterTemplateUpdater` folder, and then double-click the `ReporterTemplateUpdater.exe` file.

   Tip! As an alternative, open the Reporter Template Updater utility from the **Start** menu: click **Start** > **SCIEX** > **Reporter Template Updater**.

   The Reporter Template Updater utility opens with the **Update templates** procedure selected.
2. To select the templates to update, click **Browse**, and then click **Open**.
3. To start the update procedure, click **Update**.
4. Wait until the update procedure is completed, and then close the utility.

Remove the Component for Office 2024

CAUTION:

If Microsoft Office 2024 is installed, then we do not recommend the removal of AB SCIEX Reporter Template Suite 3.2 Component for Office 2024. If the Component for Office 2024 is removed, then the templates for the Reporter software cannot be edited on systems that have Microsoft Office 2024 installed.

For systems that have a different version of Microsoft Office installed: If the Component for Office 2024 must be removed, then restore the templates to a format that is compatible with AB SCIEX Reporter Template Suite 3.2, and then remove the Component for Office 2024.

1. Log on to the computer as a user with Administrator privileges.
2. Close the Reporter software and all of the open templates.
3. Restore the templates for the Reporter software that were updated previously to a format that is compatible with AB SCIEX Reporter Template Suite 3.2. Refer to the section: Restore the Templates for the Reporter Software.
4. Open Control Panel, and then open Programs and Features.
5. Click **AB SCIEX Reporter Template Suite 3.2 Component for Office 2024**, and then click **Uninstall**.
6. Obey the on-screen instructions.
7. Close Control Panel.

Restore the Templates for the Reporter Software

Before AB SCIEX Reporter Template Suite 3.2 Component for Office 2024 is removed, use the Reporter Template Updater utility to restore the templates for the Reporter software to a compatible format that can be edited in AB SCIEX Reporter Template Suite 3.2.

1. Go to the `C:\Program Files\SCIEX\ReporterTemplateUpdater` folder, and then double-click the `ReporterTemplateUpdater.exe` file.

   Tip! As an alternative, open the Reporter Template Updater utility from the **Start** menu: click **Start** > **SCIEX** > **Reporter Template Updater**.

   The Reporter Template Updater utility opens with the **Update templates** procedure selected.

2. In the procedure list, select **Restore templates**.

3. Click **Browse**, select the templates for the Reporter software to restore, and then click **Open**.
4. Click **Restore**.
5. Wait until the procedure is completed, and then close the utility.

Changed and Added Files

AB SCIEX Reporter Template Suite 3.2 Component for Office 2024 makes changes to these files in the `C:\Program Files\AB SCIEX\` folder:

- `ReporterOfficeAddins\TemplateContentControlManager`
  - `ABSciexReporterCommandsSchema.dll` (updated)
  - `BlankTemplate.docx` (updated)
  - `Microsoft.Office.Tools.Common.dll` (new)
  - `Microsoft.Office.Tools.Common.v4.0.Utilities.dll` (updated)
  - `Microsoft.Office.Tools.dll` (new)
  - `Microsoft.Office.Tools.v4.0.Framework.dll` (new)
  - `Microsoft.Office.Tools.Word.dll` (new)
  - `Microsoft.VisualStudio.Tools.Applications.Runtime.dll` (new)
  - `AB SCIEX Reporter Template Suite 3.2 Component for Office 2024 Release Notes.pdf` (new)
  - `TemplateContentControlManager.dll` (updated)
  - `TemplateContentControlManager.dll.manoifest` (updated)
  - `TemplateContentControlManager.vsto` (updated)
- `ReporterTemplateUpdater` (folder and all of the files are new)
  - `DocumentFormat.OpenXml.dll`
  - `ReporterTemplateUpdater.exe`
  - `ReporterTemplateUpdater.exe.config`

Contact Us

Addresses



AB Sciex Pte. Ltd.
Blk33, #04-06 Marsiling Industrial Estate Road 3
Woodlands Central Industrial Estate, Singapore 739256

SCIEX Headquarters

AB Sciex LLC
250 Forest Street
Marlborough, MA 01752
USA

Customer Training

- Global: sciex.com/contact-us

Online Learning Center

- SCIEX Now Learning Hub

SCIEX Support

SCIEX and its representatives have a global staff of fully-trained service and technical specialists. They can supply answers to questions about the system or any technical issues that might occur. For more information, go to the SCIEX website at sciex.com or use one of the links that follow to contact us.

- sciex.com/contact-us
- sciex.com/request-support

Cybersecurity

For the latest guidance on cybersecurity for SCIEX products, visit sciex.com/productsecurity.

Documentation

This version of the document supersedes all of the previous versions of this document.

To find software product documentation, refer to the release notes or software installation guide that comes with the software.

To find hardware product documentation, refer to the documentation that comes with the system or component.

The latest versions of the documentation are available on the SCIEX website, at sciex.com/customer-documents.

Note: To request a free, printed version of this document, contact sciex.com/contact-us.

This document is provided to customers who have purchased SCIEX equipment to use in the operation of such SCIEX equipment. This document is copyright protected and any reproduction of this document or any part of this document is strictly prohibited, except as SCIEX may authorize in writing.

Software that may be described in this document is furnished under a license agreement. It is against the law to copy, modify, or distribute the software on any medium, except as specifically allowed in the license agreement. Furthermore, the license agreement may prohibit the software from being disassembled, reverse engineered, or decompiled for any purpose. Warranties are as stated therein.

Portions of this document may make reference to other manufacturers and/or their products, which may contain parts whose names are registered as trademarks and/or function as trademarks of their respective owners. Any such use is intended only to designate such products as those manufacturers' products and does not imply any right and/or license to use or permit others to use such manufacturers' and/or their product names as trademarks.

SCIEX warranties are limited to those express warranties provided at the time of sale or license of its products and are the sole and exclusive representations, warranties, and obligations of SCIEX. SCIEX makes no other warranty of any kind whatsoever, expressed or implied, including without limitation, warranties of merchantability or fitness for a particular purpose, whether arising from a statute or otherwise in law or from a course of dealing or usage of trade, all of which are expressly disclaimed, and assumes no responsibility or contingent liability, including indirect or consequential damages, for any use by the purchaser or for any adverse circumstances arising therefrom.

The SCIEX clinical diagnostic portfolio is For *In Vitro* Diagnostic Use. Product(s) not available in all countries. For information on availability, please contact your local representative or refer to sciex.com/diagnostics. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd., or their respective owners, in the United States and/or certain other countries (see sciex.com/trademarks).

AB Sciex™ is being used under license.

© 2026 DH Tech. Dev. Pte. Ltd.