

For research use only. Not for use in diagnostic procedures.

ekspert ultraLC 100-XL

Release: 1.0 SP#1
Sep 2012

AB SCIEX Headquarters
500 Old Connecticut Path
Framingham, MA 01701
U.S.A

Technical Support: 1 800-952-4716 or 877-740-2129
E-mail: Support@absciex.com
Fax: 1 650-627-2739
Internet: <http://www.absciex.com>

Copyright notices

ekspert ultraLC 100-XL
Copyright(c) 2012 AB SCIEX. All rights reserved.

What is in this document

Required firmware.
Redistribution and copies of the software.
Minimum PC configuration.
Scanned for viruses.
Uninstall notes.

Minimum required firmware:

ekspert 100-XL Autosampler:	Flash	0950.110-1.02 (version 1.02)
	Template 1	0013.115-01
	Template 2	0013.116-01
ekspert 100 Pump:	Flash	0760.100-1.12 (version 1.12)
ekspert 100 Column Oven:	Flash	0886.100-2.04 (version 2.04)

Redistribution and copies of the software:

You can make as many copies of the software as you wish.
You are also allowed to redistribute copies of the software.

Minimum pc configuration:

ekspert ultraLC 100-XL runs on
Microsoft Window XP-Pro (SP2) operating systems and Microsoft Windows 7 operating systems
with .Net framework 4.0 (client profile) installed. Furthermore, Analyst 1.6 must be installed.
The following PC configuration is the minimum for ekspert ultraLC 100-XL:

- Pentium III 1.0 GHz (1000 MHz).
- 512 Mb internal memory.
- Hard disk with 850 MB of free storage space (less when .Net framework is installed already)
- SVGA video driver - 1024 x 768; 16 bit "high color";
- Serial communication port for the system controlled.
- Mouse.
- CD-ROM drive.

It is strongly recommended to install the software on a local hard disk.

1024 Mb internal memory is strongly recommended for a PC that also runs other tasks like a virus scanner, remot software and/or printer tasks.

Scanned for viruses:

The software is scanned for viruses by means of

- ESET NOD32 Antivirus version 4.2.40.0
- Virus signature database: 7450 (20120905)

- Update module: 1040 (20120313)
- Antivirus and antispyware scanner module: 1365 (20120724)
- Advanced heuristics module: 1134 (20120903)
- Archive support module: 1150 (20120828)
- Cleaner module: 1057 (20120626)
- Anti-Stealth support module: 1032 (20120806)
- ESET SysInspector module: 1214 (20100127)
- Self-defense support module: 1018 (20100812)
- Real-time file system protection module: 1006 (20110921)

Uninstall notes:

When uninstalling this software, only those files are removed that were installed.